

**Manual coarse manipulator suitable for a wide range of applications.**

This manual coarse manipulator features a classical, performance-proven design. Different component combinations enable a variety of installations, ensuring wide applicability for installation of liquid hydraulic micromanipulators. Responsive manual operation and good economy have made the MMN-1 a long-term success with pioneering researchers.

Specifications:

Accessory included		Allen Wrench
Movement range	Coarse	X 30 mm, Y 30 mm, Z 30 mm
Dimensions/Weight		W 95 × D 90 × H 110 mm, 0.47 kg

BioMedical Instruments Trading

Hofheimer Straße 63
D- 65719 Hofheim
Phone: +49 6192 2 91 88 66
Fax: +49 6192 90 13 98

Email: biomedical-instruments@science-products.com



in collaboration with:



Science Products Trading SPT AG
Käferholzstr. 142
CH- 4058 Basel
Phone +41 43 488 05 61
Fax: +41 43 488 05 62