



Stereotaxic Instrument for Mice

Based on our long-established SR series, this stereotaxic instrument is designed specifically for mice and is equipped with various functions to fix the head more firmly. The auxiliary ear bar allows clamping actions to be smoothly manipulated with just one hand, letting the operator fix the external auditory canal and check the fixation firmness with the fingertips. The mouth and nose fixing clamp features a slim, compact design appropriate for the mouth of a small mouse, making it easy to confirm that the teeth are fixed firmly. In addition, the instrument allows easy back/forward and up/down position adjustment. An 18.7mm AP frame bar is provided, to which many different types of accessories can be attached besides the attached micromanipulator SM-15. The SR-5M has one AP frame bar, while the SR-6M has two.

Specifications:

Accessories	SM-15 Stereotaxic Micromanipulator, EB-5N Auxiliary Mouse Ear Bar
Dimension, Weight	SR-5M: W 400 × D 300 × H 110 mm, 8.75 kg SR-6M: W 400 × D 300 × H 110 mm, 9.2 kg

BioMedical Instruments Trading

Hofheimer Straße 63
D- 65719 Hofheim
Phone: +49 6192 2 91 88 66
Fax: +49 6192 90 13 98

Email: biomedical-instruments@science-products.com



in collaboration with:

Science Products Trading SPT AG
Käferholzstr. 142
CH- 4058 Basel
Phone +41 43 488 05 61
Fax: +41 43 488 05 62

